



Product Code	51011
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2529115
CND	C9005
EAN/UPC	8023279510119
PARAF	977817550
GTIN	08023279510119
GMDN	13755

## Description

COLOURED TRAD STETHOSCOPES - LATEX FREE

- TRAD double head - red

Coloured aluminium stethoscopes, lightweight and reliable.  
Chestpiece diaphragm Ø 43.5 mm, bell Ø 30 mm external sizes.

Brass binaural.

Tubing: length 62 cm, int. Ø 4 mm, ext. Ø 8 mm.

Supplied in cardboard box.

Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.