



|               |                |
|---------------|----------------|
| Product Code  | 51012          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2513132        |
| CND        | C9005          |
| EAN/UPC    | 8023279510126  |
| PARAF      | 977817562      |
| GTIN       | 08023279510126 |
| GMDN       | 13755          |

## Description

COLOURED TRAD STETHOSCOPES - LATEX FREE

- TRAD double head - blue

Coloured aluminium stethoscopes, lightweight and reliable.  
Chestpiece diaphragm Ø 43.5 mm, bell Ø 30 mm external sizes.

Brass binaural.

Tubing: length 62 cm, int. Ø 4 mm, ext. Ø 8 mm.

Supplied in cardboard box.

Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.