

NEONATAL REUSABLE SpO2 PROBE for 32902, 35100-1, 35151-2 - need cable



Product Code	35107
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	1811103
CND	C900301
EAN/UPC	6945040101352
UMDNS	37808
PARAF	987475959

Description

Neonatal probe 3-10 kg - reusable

Need cable:
For use with device code 35100, cable code 35108 is required
For use with device code 32902, cable code 32917 is required