

HEINE GAMMA 3.1 PULSE STETHOSCOPE - black - M-000.09.941



Product Code	51250
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 4053755136988

Description

HEINE GAMMA STETHOSCOPES - black

GAMMA 3.1 pulse - single flat chest piece

HEINE code: M-000.09.941

Weight: 110 g

The new HEINE GAMMA stethoscopes have been given a complete overhaul. Particular emphasis was placed on an excellent ergonomic design, with angled binaural and optimised ear olives. Featuring a solid aluminium chest piece, this extremely lightweight model offers excellent acoustics and proven reliability for use in practices, hospitals, care or study settings.

- single-tube and high-quality membrane for optimum sound transmission
- durable and high-quality
- optimum patient accessibility due to angled connector on the chest piece
- non-chill rims for patient comfort
- extremely comfortable to wear

Tubing length: 70.5 cm

Chest piece: Ø 48 mm

Made in Germany.