

# HEINE GAMMA 3.3 PAEDIATRIC STETHOSCOPE - black - M-000.09.943



Product Code	51256
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2788852
CND	C9005
EAN/UPC	4053755137015

## Description

HEINE GAMMA STETHOSCOPES - black

GAMMA 3.3 paediatric - double chest piece  
HEINE code: M-000.09.943  
Weight: 115 g  
The new HEINE GAMMA stethoscopes have been given a complete overhaul. Particular emphasis was placed on an excellent ergonomic design, with angled binaural and optimised ear olives. Featuring a solid aluminium chest piece, this extremely lightweight model offers excellent acoustics and proven reliability for use in practices, hospitals, care or study settings.  
- single-tube and high-quality membrane for optimum sound transmission  
- switchable double chest piece for high and low frequencies  
- open bell diameter 30 mm  
- durable and high-quality  
- optimum patient accessibility due to angled connector on the chest piece  
- non-chill rims for patient comfort  
- extremely comfortable to wear  
Tubing length: 70.5 cm  
Chest piece: Ø 34 mm  
Made in Germany.