

HEINE GAMMA 3.2 ACOUSTIC STETHOSCOPE X SERIES - black - M-000.09.947



Product Code	51258
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 4053755203048

Description

HEINE GAMMA STETHOSCOPES - black

GAMMA 3.2 acoustic X series - double chest piece
HEINE code: M-000.09.947
Weight: 120 g

The new HEINE GAMMA stethoscopes have been given a complete overhaul. Particular emphasis was placed on an excellent ergonomic design, with angled binaural and optimised ear olives. Featuring a solid aluminium chest piece, this extremely lightweight model offers excellent acoustics and proven reliability for use in practices, hospitals, care or study settings.

- single-tube and high-quality membrane for optimum sound transmission
- switchable double chest piece for high and low frequencies
- open bell diameter 30 mm
- durable and high-quality
- optimum patient accessibility due to angled connector on the chest piece
- non-chill rims for patient comfort
- extremely comfortable to wear

Tubing length: 70.5 cm
Chest piece: Ø 48 mm
Made in Germany.