

RIESTER DUPLEX ALUMINIUM STETHOSCOPE - adult - slate grey



Product Code	51301
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396118492
---------	---------------

Description

RIESTER DUPLEX® ALUMINIUM STETHOSCOPES

- Riemer code: 4001-02
- Duplex® - slate grey
- Membrane: Ø 48 mm
- Length: 77 cm
- Weight: 151 g

Riemer DUPLEX® stethoscopes combine first-class quality at an attractive price.

They offer outstanding acoustic performance.

Special membrane for perfect auscultation with particularly flat and soft non-chill rims for improved adaptation on the skin.

Bell Ø 36 mm.

Supplied with a pair of replacement eartips and a replacement membrane.

Double-headed chest-piece in lightweight aluminium, galvanized.

Made in Germany.