

# RIESTER DUPLEX ALUMINIUM STETHOSCOPE - neonatal - blue



Product Code	51310
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	4045396006348
PARAF	951304930

## Description

### RIESTER DUPLEX® ALUMINIUM STETHOSCOPES

- Riester code: 4051
- DUPLEX® NEONATAL - blue
- Membrane: Ø 28 mm
- Length: 76 cm
- Weight: 106 g

Riester DUPLEX® stethoscopes combine first-class quality at an attractive price.

They offer outstanding acoustic performance.

Special membrane for perfect auscultation with particularly flat and soft non-chill rims for improved adaptation on the skin.

Bell Ø 36 mm.

Supplied with a pair of replacement eartips and a replacement membrane.

Double-headed chest-piece in lightweight aluminium, galvanized.

Made in Germany.