

RIESTER CARDIOPHON 2.0 STAINLESS STEEL STETHOSCOPE - black - 4240-01



Product Code	51314
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2934754
CND	C9005
EAN/UPC	4045396195547

Description

RIESTER CARDIOPHON 2.0 STETHOSCOPES

- Cardiophon 2.0 - black
- Riester code: 4240-01
- Membrane mm: Ø 32/44
- Length cm: 71
- Weight g: 204

Developed for cardiology, setting a new standard for precision and acoustic performance. The completely reconstructed chest-piece, sound transmission and binaural, together with several additional innovations, deliver a superior stethoscope excelling in the frequency range between 200 Hz and 500 Hz, which is particularly important for auscultation of heart sounds. Suitable for adults and children.

- innovative acoustic system with extraordinary auscultation for all frequency ranges
- high-precision double chest-piece made of stainless steel
- bilateral special membrane for precise acoustics, Ø 44 and 32 mm with novel non-chill rims for improved skin comfort
- the binaural is anatomically shaped and equipped with internal multiple springs for adjusting contact pressure
- extra-soft, replaceable and pivotable ear-tips for perfect sealing of the ear canal and enhanced wear comfort
- Y-tube with separate sound-conducting pathways for each ear, with excellent audio-transmission, which allows to hear the subtle characteristics of heart sounds and murmurs more distinctly.

Latex-free. Supplied with 2 pairs of ear-tips, one replacement membrane for each and name tag.

Made in Germany.