

RIESTER CARDIOPHON 2.0 STAINLESS STEEL STETHOSCOPE - blue - 4240-03



| | |
|---------------|----------------|
| Product Code | 51316 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2934754 |
| CND | C9005 |
| EAN/UPC | 4045396195561 |

Description

RIESTER CARDIOPHON 2.0 STETHOSCOPES

- Cardiophon 2.0 - blue
- Riemer code: 4240-03
- Membrane mm: Ø 32/44
- Length cm: 71
- Weight g: 204

Developed for cardiology, setting a new standard for precision and acoustic performance. The completely reconstructed chest-piece, sound transmission and binaural, together with several additional innovations, deliver a superior stethoscope excelling in the frequency range between 200 Hz and 500 Hz, which is particularly important for auscultation of heart sounds. Suitable for adults and children.

- innovative acoustic system with extraordinary auscultation for all frequency ranges
- high-precision double chest-piece made of stainless steel
- bilateral special membrane for precise acoustics, Ø 44 and 32 mm with novel non-chill rims for improved skin comfort
- the binaural is anatomically shaped and equipped with internal multiple springs for adjusting contact pressure
- extra-soft, replaceable and pivotable ear-tips for perfect sealing of the ear canal and enhanced wear comfort
- Y-tube with separate sound-conducting pathways for each ear, with excellent audio-transmission, which allows to hear the subtle characteristics of heart sounds and murmurs more distinctly.

Latex-free. Supplied with 2 pairs of ear-tips, one replacement membrane for each and name tag.

Made in Germany.