

# RIESTER RI-SONIC ELECTRONIC STETHOSCOPE - USB-A plug - 4300



Product Code	51321
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2275230
CND	V0399
EAN/UPC	4045396223523

## Description

### RI-SONIC PCP-USB ELECTRONIC STETHO- USB-A plug - 4300

The Riester Ri-sonic are electronic stethoscopes designed to provide real-time and remote auscultation.

PCP-1 plugs directly into microphone port of any smartphone, PC or tablet and is only suitable for eMurmur App.

PCP-USB stethoscopes connect directly to any smartphone, PC or tablet and record sounds via eMurmur App but they're also ideal for Telemedicine projects thanks to the integrated embedded amplifier.

- patented piezo technology for picking up body sounds with minimal background artefacts
- patients can use the stethoscope on themselves
- easy to use with a guidance of a clinician.

Ri-sonic Electronic Stethoscopes enable to listen to internal body sounds, such as heart, lung or bowel sounds, of a patient.

## Technical Specifications

Mode of operation: continuous operation

Electrical protection: Class II protection against electrical shock

Plug: USB plug

Voltage: Input: 2 V - 5 V DC

Classification: Class IIa Type BF applied part

Weight: 180 g

eMURMUM App

Advanced digital auscultation platform.

It makes possible for professionals to screen, monitor, diagnose and consult heart, lung and bowel health in clinics, hospitals and home.