

3M SINGLE-PATIENT PAEDIATRIC STETHOSCOPE - SPS-YP1010 - yellow



Product Code	51400
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 0707387784195

Description

3M™ SINGLE-PATIENT STETHOSCOPES - length 72.4 cm

- 3M code: SPS-YP1010
- Paediatric - head Ø 36.2 mm
- Box of 10

3M™ has developed a new stethoscope, which combines advanced acoustic engineering and infection prevention. It is intended for single-patient use, adult or child in an isolation area, and to be disposed when no longer needed. It helps to reduce potential risks of patient-to-patient crosscontamination, that could occur using reusable stethoscope.

Features:

- high-quality acoustics with clear and reliable sound
- soft and flexible earplugs isolate from noises
- pressure sensitive diaphragm allows to hear high and low frequency sounds
- comfortable flexible headset and angled ear tubes provide an excellent seal
- flexible tubing, to easily move and position chestpiece
- easy-grip moulded chestpiece for multiple hand positions