

3M SINGLE-PATIENT ADULT STETHOSCOPE - SPS-YA1100 - yellow



Product Code	51415
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

Description

3M™ SINGLE-PATIENT STETHOSCOPES - length 72.4 cm

- 3M code: SPS-YA1100
- Adult - head Ø 45.2 mm
- Box of 100

3M™ has developed a new stethoscope, which combines advanced acoustic engineering and infection prevention. It is intended for single-patient use, adult or child in an isolation area, and to be disposed when no longer needed. It helps to reduce potential risks of patient-to-patient crosscontamination, that could occur using reusable stethoscope.

Features:

- high-quality acoustics with clear and reliable sound
- soft and flexible earplugs isolate from noises
- pressure sensitive diaphragm allows to hear high and low frequency sounds
- comfortable flexible headset and angled ear tubes provide an excellent seal
- flexible tubing, to easily move and position chestpiece
- easy-grip moulded chestpiece for multiple hand positions