

LITTMANN KIT 40006: DIAPHRAGM+RIM+BELL SLEEVE +EARTIPS for Classic II SE - grey - blister



Product Code	51465
Unit of sale	1 kit
Minimum order	1
Type	Spare part
Class	

GTIN 00707387772970

Description

Kit accessories 3M™ Littmann®

TUNABLE DIAPHRAGM AND RIM + NON CHILL BELL SLEEVES + EARTIPS

- 3M Code: 40006
- For: Classic II S.E.
- Colour: grey
- blister