

LITTMANN KIT 40021: DIAPHRAGM+RIM+BELL SLEEVE +EARTIPS for Lightweight - light brown - blister



Product Code	51467
--------------	-------

Unit of sale	1 kit
--------------	-------

Minimum order	1
---------------	---

Type	Spare part
------	------------

Class	
-------	--

GTIN	00707387773069
------	----------------

Description

Kit accessories 3M™ Littmann®

TUNABLE DIAPHRAGM AND RIM + NON CHILL BELL SLEEVES + EARTIPS

- 3M Code: 40021
- For: Lightweight
- Colour: light brown
- blister