

# LITTMANN KIT 40016: 2 DIAPHRAGMS+RIM+EARTIPS for Classic III, Cardiology IV - black - blister



Product Code 51482

Unit of sale 1 kit

Minimum order 1

Type Spare part

Class

GTIN 00707387783969

## Description

Kit accessories 3M™ Littmann®

2 TUNABLE DIAPHRAGMS + RIM ASSEMBLY + EARTIPS

- 3M Code: 40016
- For: Classic III™, Cardiology IV™, Core
- Colour: black
- blister