

LITTMANN KIT 40017: 2 DIAPHRAGMS+RIM+EARTIPS for Classic III, Cardiology IV - grey - blister



Product Code 51483

Unit of sale 1 kit

Minimum order 1

Type Spare part

Class

GTIN 00707387784010

Description

Kit accessories 3M™ Littmann®

2 TUNABLE DIAPHRAGMS + RIM ASSEMBLY + EARTIPS

- 3M Code: 40017
- For: Classic III™, Cardiology IV™, Core
- Colour: grey
- blister