

# IMPERMEABLE, ABSORBENT ROLL 90 cm x 50 m - 3500 ml/m<sup>2</sup>



Product Code	25167
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2184860
CND	T030104
EAN/UPC	8023279251678

## Description

Absorbent rolls. Can be cut at desired size.

- Size: 90 cm x 50 m
- Absorbency: 3,500 ml/m<sup>2</sup>

Roll mats to be positioned on the floor of the operating room or near the sink in the pre-operating theatre for holding liquids that could fall on the floor.

Mats have high absorbent capacity and optimal speed absorbency.

The mats distribute the absorbed liquid by capillary method and maintain it at his internal structure avoiding the lifts on surface.

The special film at rear side blocks the mat on the floor and keeps it stopped without using adhesive matters, which could release glue traces on the floor.

They are perfectly adherent to the ground, avoiding risk of accidents.