

CHECKME LITE HEALTH MONITOR - WIRELESS



Product Code	35185
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2056594
CND	Z1203020202
EAN/UPC	6934440100041
PARAF	977815582

Description

CHECKME™ LITE HEALTH MONITOR - WIRELESS

Innovative all-in-one device with multiple functionalities in one compact and portable unit. It allows simultaneous measurement of ECG, heart rate, QRS duration, SpO2, pulse rate and pulse index. With colour touch screen, well-organized user interface and audio instructions, it is very intuitive and easy to use. Clinically approved accuracy.

- supports both cable and cable-free measurement of ECG and SpO2 (optional cable)
- checks, captures and identifies abnormal ECG in 20 sec.
- professional ECG algorithm provides heart rate, QRS duration, ST segment and rhythm analysis
- checks blood oxygenation with built-in SpO2 sensor in 20 sec
- high-quality ECG waveform can be saved and replayed
- professional PC data management software for review and generating reports (compatible with Windows OS)
- ViHealth App.

Manual: GB, IT, ES, DE, FR
Audio instructions: GB, PL, DE, IT, FR, ES, HU
PC Software: GB

ECG + HR + SpO2 + QRS + Pulse Rate + PI
Size: 88x56x13 mm
Weight: 64 g

Standard accessories

PC software: GB
ECG cable
ECG electrodes
Neck strap
Charging cable, micro D