

FLAVIA LIFT ARMCHAIR 2 motors - burgundy



Product Code	45188
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279451887
PARAF	977817028

Description

FLAVIA LIFT ARMCHAIRS WITH 2 MOTORS
2 motors lift armchair with standard or extra lifting footrest: it can reach the bed position providing a relaxing experience. Its backrest and footrest can operate separately to any position. Extra lifting footrest allows 8 cm lifting adjustment higher than flat position.

- Color: burgundy
- 2 motors
- Fireproof: NO
- Footrest: standard

Lift armchairs that feature a powered lifting mechanism that pushes the chair up from its base to assist the user to a standing position.

Power supply: 100-240 V - 50/60 Hz.

Multilingual manual: GB, FR, IT, ES.

Technical Specifications

- Seat
Width: 47 cm
Depth: 50 cm
Height: 47 cm
- Backrest height: 74 cm
- Chair
Width: 78 cm
Lenght (reclined): 171.5 cm
- Load: 150 Kg (331 lbs)