

STREP-A TEST - dual colour - strip



Product Code	24519
Unit of sale	box of 50 pcs
Minimum order	1
Type	IVD - Professional Use
Class	
RDM (NSIS)	1702284
CND	W0105090103
EAN/UPC	8023279245196
PARAF	973922469

Description

STREP-A DUAL COLOUR TEST - strip - box of 50 tests

Dual colour rapid test for Strep A uses the antibody labelled colour particles coated at two separate locations in the test device. Optical immunoassay technology allows greater sensitivity and less false negative results without affecting specificity.

Standard procedure and easy to view results. If Group A Streptococcus is present in the sample, a visible red Test Line will appear to indicate a positive result and a blue Control Line will also appear to indicate the test is valid.

Kit includes 50 test strips, 50 tubes, 50 swabs, 2 reagents and 2 control solutions.

WARNING: SHORT EXPIRY DATE FOR TESTS. Most tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

Technical Specifications

Relative sensitivity: 95.1% (95%CI: 89.6-98.2%)
Relative specificity: 97.8% (95%CI: 95.8-99%)
Overall accuracy: 97.1% (95%CI: 95.3-98.4%)