

ARMFIT+ BLOOD PRESSURE MONITOR WITH ECG



Product Code	35198
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2555883
CND	Z1203020501
EAN/UPC	6934440102212
UMDNS	16173
PARAF	981509375

Description

ARMFIT+ BLOOD PRESSURE MONITOR WITH ECG

Armfit+ is intended to be used for real time monitoring, measuring, displaying, reviewing and storing ECG and blood pressure variations, at home or healthcare facilities.

ECG can be recorded without cable only.

Built-in Bluetooth wireless connectivity allows to share measurement records to mobile phones and tablets.

Supplied with one adult cuff, charging cable, quick guide and multilingual manual in GB, FR, IT, ES, DE.

Technical Specifications

Blood Pressure Measurements

Technology: Oscillometric Method

Pressure measurement range: 0-300 mmHg

Accuracy: ± 3 mmHg

Pulse rate range: 40 to 200 /min

Pulse rate accuracy: ± 2 /min

Clinical accuracy: Meet IEC80601-2-30

Blood pressure records: 50

ECG Recording

Lead type: Integrated ECG electrodes

Lead set: Lead I, Lead II, Chest Lead

ECG measurement: 30s recording

Heart rate range: 30-250/min ± 2 /min

ECG records: 10

Device

Size - weight (main unit): 135×45×20 mm - 240 g

Cuff size: Adult cuff: 22-42cm

Wireless: Built-in Bluetooth 4.0

Charge input: Micro USB, DC5V

Battery type: Rechargeable lithium-polymer

Battery run time: 500 measurements

Charge time: 2 hours

Classifications

Protection: Type BF

Dust & water resistance: IP22

