

# S/S STRAIGHT SCISSORS with belt clip - blunt/sharp - 14 cm



Product Code	20520
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2775713
CND	L01040199
EAN/UPC	8023279205206
GTIN	08023279205206
GMDN	13481

## Description

Stainless steel nurse straight scissor - 14 cm

- Model: Standard with belt clip

Powder coated instruments are not autoclavable, can be sterilized in water at maximum 110°C