

FOLDED NAPKINS - 33x45 cm blue



|               |                |
|---------------|----------------|
| Product Code  | 25226          |
| Unit of sale  | box of 500 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2442610       |
| CND        | T0399         |
| EAN/UPC    | 8023279252262 |
| PARAF      | 987472560     |

Description

FOLDED NAPKINS - pack of 500 pcs  
Napkins 1+1 (1 ply cellulose + 1 ply polyethylene)

- Napkins 33x45 cm - blue

Pre-coated folded napkins, made of 1 ply of pure cellulose wadding laminated with 1 ply of polyethylene.  
High absorbency and waterproofness to assure high comfort.

Made in Italy.