

SpO2 ADULT "Y" EAR PROBE for DATEX-OHMEDA - 3.0 m cable



Product Code	35235
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100542

Description

OEM - Brand: Datex-Ohmeda
N: 3
Sensor type: Adult Y-type ear sensor
Cable Length: 3.0 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Datex-Ohmeda: AS/3 ,CS/3, CH/S, CH/2S, CH-RS, SAF-F, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satlite Trans, Satlile, Trans Plus, System 5, ULTS00