

SpO2 ADULT PROBE for GE DATEX-OHMEDA - 3.0 m cable



Product Code	35238
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100337

Description

OEM - Brand: GE Datex-Ohmeda
N: 1
Sensor type: Adult Finger
Cable Length: 3.0 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: GE (Datex-Ohmeda): S/5, i4 Modular Monitors with Eseries modules, Corometrics® 259 with TruSignal Ohmeda Technology