



Product Code	25271
Unit of sale	box of 2 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2192129
CND	H90020199
PARAF	930048451
GTIN	07610182013076
GMDN	32207

Description

3M™ STERI-STRIP™ - STERILE

- 3M™ STERI-STRIP™ - 6x75 mm - white
- 3M™ code: 1541NP

3M™ Steri-Strip™ reinforced adhesive skin closures are porous non-woven adhesive strips indicated for use as a skin closure device in the treatment of lacerations and surgical incisions.

Designed to let the skin breathe and reinforced with filaments for strength.

- fast and easy to apply
 - comfortable to wear compared to sutures
 - non-invasive and sterile design helps to reduce risk of infections
- Pressure sensitive hypoallergenic adhesive.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, GR.

Box: FR, IT, DE.