

# SpO2 ADULT/NEONATAL WRAP PROBE for PHILIPS - 1.6 m cable



|               |                |
|---------------|----------------|
| Product Code  | 35276          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II B           |

|            |               |
|------------|---------------|
| CH-REP     | No            |
| RDM (NSIS) | 2621198       |
| CND        | C900301       |
| EAN/UPC    | 6900234100511 |

## Description

OEM - Brand: Philips  
N: 4  
Sensor type: Adult/neonatal wrap  
Cable Lenght: 1.6 m  
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Philips: M3, M4, MP20, MP30, MP40, MP50, MP70, MP80, VM4, VM6, VM8