

SpO2 ADULT SOFT PROBE for PHILIPS - 3.0 m cable



Product Code	35283
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100474

Description

OEM - Brand: Philips
N: 2
Sensor type: Adult soft-finger
Cable Lenght: 3.0 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Philips: 1205, M1205A, 78352C, 78834C, 78354C, M1020A, M1025A, M1025B, M1350B, M1772A/B, M1732A/B, M2475B