

SpO2 ADULT "Y" EAR PROBE for SIEMENS/DRAGER - 1.6 m cable



Product Code	35287
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100610

Description

OEM - Brand: Siemens/Drager

N: 3

Sensor type: Adult Y-type ear

Cable Length: 1.6 m

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Siemens/Drager: SC6002XL, 7000, SC8000, SC9000XL