

DAE-AED DOUBLE-SIDED ALUMINIUM SIGN 34x36 cm
for defibrillators



Product Code	35327
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8034028012211
---------	---------------

Description

DAE-AED DOUBLE-SIDED ALUMINIUM SIGN - 34x36 cm

Signs in aluminium foil with thickness 5/10.
Prepainted white with silkscreen in green colour.
Pre-drilled for wall mounting.

Weight: 350 g