

CABINET FOR DEFIBRILLATORS - outdoor use



Product Code	35334
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8034028012174

Description

METAL CABINET FOR DEFIBRILLATORS - OUTDOOR USE - yellow

Not suitable for marine or corrosive environments. Special coated cabinets resistant to mentioned environments available on request in 30 days.

Water tight welding on the sides, water and dust proof gaskets.
Made of carbon steel, thickness 15/10.
Adjustable internal shelf.
Door with transparent acrylic window.

External size: 43x21xh 48 cm
Internal size: 42x20xh 45 cm