

SANIBEL ADI FLANGED INFANT EAR TIP 5-8 mm - yellow



Product Code	53517
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 8023279535174

Description

SANIBEL™ ADI FLANGED INFANT EAR TIPS

Ear tip 5-8 mm - yellow

SANIBEL™ ear tips were designed for use with acoustic immittance instruments and otoacoustic emission system. Made of a soft and sticky biocompatible material, these ear tips greatly improve the ease of obtaining and keeping a good seal, allowing the tester to quickly move on to other tasks. Available in different sizes, with different shapes (mushroom and flanged ear tips for infants). The complete range makes it possible to accommodate all ear canal shapes and sizes. All Sanibel™ ear tips are latex free.

Technical Specifications

VERIFIED COMPATIBILITY FOR SANIBEL™ ADI

- Amplivox® Otowave™ with clear (Type 1) probe
- Interacoustics® AT235™ & AT235h™ with clear probe tip adapter, AA222™ & AA220™ with clear probe tip adapter, Eclipse™ and MT10™ with clear probe tip, Titan™
- Maico® MI24™, MI26™, MI34™, MI44™ with clear probe tip adapter, easyTimp™ - Otowave™ - Ero-scan plus™ - easyScreen™
- MEDRx® Otowave™
- Grason-Stadler®: Corti™