

# MAICO MA27 SCREENING AUDIOMETER - air conduction



Product Code	53549
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2233028
CND	Z121401
EAN/UPC	4260176120865

## Description

### MAICO MA27 SCREENING AUDIOMETER

Air conduction audiometer for basic screening with rigid and lightweight case. It has 11 frequencies, built-in carrying handle and compartment to store headphones and accessories. Tough, accurate and reliable with large fingertip control knobs for easy frequency and hearing levels adjustments.

Selectable display languages: GB, FR, IT, ES, DE.  
Supplied with AC headset DD45, audiogram pad and power supply.

Multilingual manual: GB, FR, IT, ES, DE, PL, RU.

## Technical Specifications

Display: 38.1x76.2 mm, monochrome  
Level steps: 5 or 1 dB  
Hearing level range: -10 dB HL to 100 dB HL  
Test signal: pure tone and warble tone (pulsed and continuous)  
Test frequencies: 125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000 Hz  
Size - Weight: 25.5x37x15 cm - 2.4 kg  
Power supply: 100-240VAC  $\pm$  10 %, 50/60 Hz, 0.6 A  
Standards: IEC 60601-1:2012, AAMI ES60601-1:2005+A2+A1, CAN/CSA-C22.2 No. 60601-1:14 Type B IEC 60601-1-2:2014, IEC 60645-1:2017 Type 4, ANSI/ASA S3.6-2018 Type 4  
Classification: Class IIa (EU 2017/745), Medical Device