

AMPLIVOX 240 DIAGNOSTIC AUDIOMETER - air, bone, mask conduction



Product Code	53551
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2551529
CND	Z121401
EAN/UPC	8023279535518

Description

AMPLIVOX 240 DIAGNOSTIC AUDIOMETER - air, bone, masking

Cost-effective diagnostic audiometer with USB connectivity, ideal for desktop and mobile use.

- 11 test frequencies: 125 to 8,000 Hz (AC), 250 to 8,000 Hz (BC)
- -10 to 120 dBHL range (AC), -10 to 70 dBHL (BC), 90 dBHL max (250 to 4,000 Hz) insert masking output
- threshold retention function
- single, pulsed or warble tone preset
- automatic and manual AC testing
- optional battery power function

Supplied with audiometric headset, bone conductor headset, patient response switch, NOAH and ampliSuite software, 50 audiogram cards, main power adaptor, calibration certificate, carrying case and manual GB, IT, ES. On request FR, PT, PL, RO, FI.

Weight: 750 g

Dimension: 270x165x60 mm

Operating power: 100-240 V, 50-60 Hz