



Product Code	53564
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8054187380747
---------	---------------

Description

SCREENING PURE TONE AUDIOMETRY (PTA) LICENSE WITH PRE-CALIBRATED HEADSET

To conduct pure tone air conduction audiometry.

Technical Specifications

- Audiometry
- Intensity: -10 to 100 db HL
- Freq. range: 125-8000 Hz
- Input: pure tone
- Output: AC