

PHARMADOCT TWEEZER - carton of 12 boxes



Product Code	25363
Unit of sale	1 carton
Minimum order	1
Type	No medical device
Class	

UK-REP	Yes
CH-REP	Yes
EAN/UPC	8017990160555
PARAF	976327243

Description

PHARMADOCT: first aid, medication and wellness.

- Tweezer
- Pieces per box: 1
- Minimum order: 12 box

Graphic of Pharmadoc line boxes are subject to changes.

All products are in blister with 7 languages package (GB, FR, IT, ES, DE, PT, GR).