



Product Code	53701
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2715820
CND	Z12030102
EAN/UPC	6922166500054
PARAF	987471796

## Description

### BREATHCARE PAP DEVICE - CPAP/AUTO CPAP

The breathcare PAP devices take care of patients every night.

PAP (Positive Airway Pressure) device works as an air compressor, it constantly blows pressurized air to patient's airway, to prevent airway collapse.

The PAP devices contain CPAP and Auto CPAP.

CPAP: Continuous Positive Airway Pressure device

Auto CPAP: Automatic Continuous Positive Airway Pressure device

Immediately improve sleep quality

No trauma, safe and comfortable

User friendly with little discomfort

Quiet, noise = 28 dB

Ultra quiet thanks to 4 methods for reducing noise level.

Stable 20,000 hours life time.

Comfortable based on years of clinical data it effectively reduces the average treatment pressure and provides a comfortable therapy experience.

Connection with oxygen concentrator to supply oxygen with mask which can help ease symptoms of hypercapnia.

Internal software: GB, IT. Supplied with carrying bag, 2 filters, air hose, heated humidifier, breathing nose mask, SD card, power adapter, power cable and manual (GB, IT).

## Technical Specifications

Treatment mode: CPAP , Auto CPAP

Pressure range: 4~20 cm H<sub>2</sub>O ±0.5 cm H<sub>2</sub>O

Time-lapse boosting: 0~45 min

Alarm protection: face mask fell alarm

Operating power: 100-240 V, 50-60 Hz DC 24 V

Automatic altimetric compensation: 3500m

USM air circuit system: ultra silence air circuit system

PCM: noise removal technology

3A self adapting frequency adjustment technique

Size - weight: 27x13.5x10 cm - 1.5 kg