



Product Code	53705
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2810597
CND	R0301010201
EAN/UPC	6922166510527
UMDNS	12447
PARAF	981516228

Description

Breathing nose mask - medium

Airtight, light and soft masks adjust automatically according to different face shapes. Honeycomb vents designed to be gentle and quiet. Quick-release connector helps separate the tubing from the mask easily. Headgear is made of soft and skinfriendly fabric to ensure maximum comfort. Masks are suitable for Yuwell, Resmed, Philips, Fisher & Paykel and most of CPAP machines in the market. Therapy pressure range: 4-30 cm H₂O.
Tube connection: 22 mm

Weight: 64 g