



Product Code	53900
Unit of sale	box of 300 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	930556
CND	C020501
GMDN	35035

Description

AMBU® WhiteSensor

WhiteSensor electrodes combine versatility with quality and functionality.

The wide range of sizes, shapes, gel types and backing materials makes it easy to choose the optimal device for every situation.

Carefully selected features to match each individual situation, covering all patient needs in monitoring.

Flexible and reliable materials and adhesives, able to follow skin movements without lifting off, allow both patient comfort and signal quality.

Ambu® WhiteSensor 4540

Highly conductive wet gel system which guarantees a fast and accurate trace during 12-lead ECG procedures.

Thanks to the strong high-tack adhesive foam, time wasted due to electrode drop-off will be reduced.

The radiolucent 4 mm slot connector allows the patient to undergo X-ray and scanner procedures.

5 pieces per card.

Technical Specifications

- AMBU code: 4540
- Description: WhiteSensor 4540
- Application: Resting
- Gel type: Wet
- Backing material: Foam
- Size(mm): 32x44
- Shape: Rectangular
- Shell life* (month): 24
- Radiolucent: YES
- Connector type: 4 mm
- Minimum order: Box of 300 pcs.
- Packaging: 5 pouches of 60 pcs.