

# AMBU WHITESENSOR WS-00-S/50 ECG ELECTRODES

- snap



Product Code	53904
Unit of sale	box of 600 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	556932
CND	C020501
GTIN	28906035282694
GMDN	35035

## Description

### AMBU® WhiteSensor

WhiteSensor electrodes combine versatility with quality and functionality.

The wide range of sizes, shapes, gel types and backing materials makes it easy to choose the optimal device for every situation.

Carefully selected features to match each individual situation, covering all patient needs in monitoring.

Flexible and reliable materials and adhesives, able to follow skin movements without lifting off, allow both patient comfort and signal quality.

### Ambu® WhiteSensor WS-00-S/50

Highly conductive solid gel with strong adhesion to ensure a good signal quality during short-term ECG monitoring applications.

Thanks to the flexible foam backing material, the electrode ensures ease of use and comfort during application.

## Technical Specifications

- AMBU code: WS-00-S/50
- Description: WhiteSensor WS-00-S/50
- Application: Short-term
- Gel type: Solid
- Backing material: Foam
- Size(mm): 36x40
- Shape: Teardrop
- Shell life\* (month): 30
- Radiolucent: NO
- Connector type: Snap
- Minimum order: Box of 600 pcs.
- Packaging: 12 pouches of 50 pcs.