



Product Code	53906
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	930565
CND	C020501
GTIN	18906035283342
GMDN	35035

Description

AMBU® WhiteSensor

WhiteSensor electrodes combine versatility with quality and functionality.

The wide range of sizes, shapes, gel types and backing materials makes it easy to choose the optimal device for every situation.

Carefully selected features to match each individual situation, covering all patient needs in monitoring.

Flexible and reliable materials and adhesives, able to follow skin movements without lifting off, allow both patient comfort and signal quality.

Ambu® WhiteSensor 0315M

Small size tab electrode with solid gel for 12-lead ECG.

The fish tail tab design prevents electrode lift-offs even if the adapter kinks.

10 pieces per card.

Technical Specifications

- AMBU code: 0315M
- Description: WhiteSensor 0315M
- Application: Short-term
- Gel type: Solid
- Backing material: Polymer
- Size(mm): 35x21
- Shape: Tab
- Shell life* (month): 24
- Radiolucent: YES
- Connector type: Tab
- Minimum order: Box of 500 pcs.
- Packaging: 10 pouches of 100 pcs.