

AMBU WHITESENSOR 4500M-H ECG ELECTRODES - snap



Product Code	53910
Unit of sale	box of 300 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	993287
CND	C020501
GTIN	18906035282918
GMDN	35035

Description

AMBU® WhiteSensor

WhiteSensor electrodes combine versatility with quality and functionality.

The wide range of sizes, shapes, gel types and backing materials makes it easy to choose the optimal device for every situation.

Carefully selected features to match each individual situation, covering all patient needs in monitoring.

Flexible and reliable materials and adhesives, able to follow skin movements without lifting off, allow both patient comfort and signal quality.

Ambu® WhiteSensor 4500M-H

It features a wet gel with good adhesion to ensure a good signal quality during shortterm ECG monitoring applications.

Thanks to the flexible foam backing material, the electrode ensures ease of use and comfort during application.

Technical Specifications

- AMBU code: 4500M-H
- Description: WhiteSensor 4500M-H
- Application: Stress-test Short-term
- Gel type: Wet
- Backing material: Foam
- Size(mm): 50x48
- Shape: Round
- Shell life* (month): 24
- Radiolucent: NO
- Connector type: Snap
- Minimum order: Box of 300 pcs.
- Packaging: 5 pouches of 60 pcs.