

NYTRILE GLOVES - LONG - 300 mm - medium



Product Code	25397
Unit of sale	box of 100 pcs
Minimum order	1
Type	PPE
Class	III

UK-REP	Not required
EAN/UPC	8016629006486

Description

Nitrile gloves - medium - regular blue - box of 100

Disposable nitrile examination gloves powder-free and latexfree.
The 300 mm length and 5.7 g weight ensure a high level of protection.
The micro textured surface ensures a very good grip.
Suitable for use in hospital.

Medical device. Class I - EN 455.
PPE (Personal Protection Equipment) Cat III.
Non sterile.
EN 420, EN 374, EN 388 AQL<1.5

Multilingual box: GB, FR, IT, ES, PT, DE.