

# NYTRILE GLOVES - LONG - 300 mm - large



Product Code	25398
Unit of sale	box of 100 pcs
Minimum order	1
Type	PPE
Class	III

UK-REP	Not required
EAN/UPC	8016629006493

## Description

Nitrile gloves - large - regular blue - box of 100

Disposable nitrile examination gloves powder-free and latexfree.  
The 300 mm length and 5.7 g weight ensure a high level of protection.  
The micro textured surface ensures a very good grip.  
Suitable for use in hospital.

Medical device. Class I - EN 455.  
PPE (Personal Protection Equipment) Cat III.  
Non sterile.  
EN 420, EN 374, EN 388 AQL<1.5

Multilingual box: GB, FR, IT, ES, PT, DE.