


# ECG PATIENT CABLE 2.2 m - banana - compatible Fukuda Denshi



ECG connector  
Connettore ECG  
Connecteur ECG

Banana

Product Code	54008
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2341763
CND	Z12050380
EAN/UPC	8023279540086

## Description

### ECG PATIENT CABLE 10 LEADS

- Leads: Banana
- Cable length (Length from connector to yoke): 2.2 m

### COMPATIBILITY (BRAND AND MODELS)

- FUKUDA DENSHI: Autocardiner & Cardimax. FCP11, FCP2101A, FCP2155, FCP2201A, FCP4101A, FX101, MCP302, FX2111, FCP2155, FCP7101, FX3010. FX328UY, FX4010, FX7101, FX7102, FX7202, FX7302, FX7402, FX7542, FX8222, FX8322, FX8328, FX328U

Compatible patient cables for electrocardiography and exercise stress tests.

In order to ensure the patient's safety under any circumstance, we have supplied all our cables with an independent defibrillation protection.

Quality is ensured by the design, the manufacturing process and the final tests performed through computerised workstations.

Standard cable length 2.2 m (from connector to yoke), chest leads 0.76 m, limb leads 1.2 m.

Few cables with 3.5 m length to improve usability in exercise stress tests.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, RO.