


ECG PATIENT CABLE 2.2 m - banana - compatible Esaote, Shiller



ECG connector
Connettore ECG
Connecteur ECG

Banana

Product Code	54011
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2341763
CND	Z12050380
EAN/UPC	8023279540116

Description

ECG PATIENT CABLE 10 LEADS

- Leads: Banana
- Cable length (Length from connector to yoke): 2.2 m

COMPATIBILITY (BRAND AND MODELS)

- ESAOTE: P80, P80six, P80power, P8000, P8000power, P8000 tele, MyCardiopad, MyCardiopad XL, AT-10 plus, AT-102 plus, Cardiovit M2015 ecc
- SCHILLER: AT-1, 1G2, -2, -2plus, -2light, -4, -101, -104, -110, -5, -6, -10, -10+, -10plus, -102, -102G2, -60, -100, -104PC, CS-6, -6/12, -12, -100, series -200 MS-2007, -2010, -2015

Compatible patient cables for electrocardiography and exercise stress tests.

In order to ensure the patient's safety under any circumstance, we have supplied all our cables with an independent defibrillation protection.

Quality is ensured by the design, the manufacturing process and the final tests performed through computerised workstations.

Standard cable length 2.2 m (from connector to yoke), chest leads 0.76 m, limb leads 1.2 m.

Few cables with 3.5 m length to improve usability in exercise stress tests.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, RO.