


ECG PATIENT CABLE 3.5 m - snap - compatible Camina, Colson, ST, others

Product Code	54020
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I



ECG connector
Connettore ECG
Connecteur ECG

Snap

RDM (NSIS)	2341763
CND	Z12050380
EAN/UPC	8023279540208

Description

ECG PATIENT CABLE 10 LEADS

- Leads: Snap
- Cable length (Length from connector to yoke): 3.5 m

COMPATIBILITY (BRAND AND MODELS)

- CAMINA-EUROCAMINA: Cardiorapid and CardioSwift series
- CARDIETTE-ELETTRONICA TRENINA: AR600, 600adv, 600view; AR1200, 1200adv, 1200view; AR2100, 2100adv, 2100view. Start 100P, 100H, 200, 200HV. Daedalus View, View HES, View Ergo. DOC 101,103; micro, microruler and microtel series; Autoruler series; Excel series: Daedalus 346
- CARDIOLINE: ECG100 and ECG200 series. All models AR600, 1200, 2100; E1, E1R, E3
- COLSON - DUPONT: Cardi-3, E1, E3 Cardiopraticque, Daedalus 346. All models of the series: AR, Cardipocket, Cardio, Excel, Autoruler, Microruler. All models with Dsub connector 15 pins
- FARUM: All models with Dsub connector 15 pins
- MEGOS: 3/12b, 3/12i; Comp1x12, Comp3x12, Comp6x12, Comp3, Comp12, Comp60 and models with Dsub connector 15 pins
- ST ELECTROMEDICINA: 531, 531i, 563, 563i. ar600adv, ar600view, ar1200adv, ar1200view

Compatible patient cables for electrocardiography and exercise stress tests.

In order to ensure the patient's safety under any circumstance, we have supplied all our cables with an independent defibrillation protection.

Quality is ensured by the design, the manufacturing process and the final tests performed through computerised workstations.

Standard cable length 2.2 m (from connector to yoke), chest leads 0.76 m, limb leads 1.2 m.

Few cables with 3.5 m length to improve usability in exercise stress tests.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, RO.