



Product Code	54100
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
EAN/UPC	4015672114012

Description

KARDIAMOBILE 1 LEAD ECG - KM-1L-E

KardiaMobile are pocket-sized personal ECGs with finger electrodes that work together with smartphone or tablet through Kardia App.

Visit the following link to view the full list of countries where the product and app are currently available:
<https://alivecor.zendesk.com/hc/en-us/articles/115015719967-Which-Kardia-products-are-available-in-my-Country>

The KardiaMobile System (device + Kardia App) displays ECG rhythms and output of ECG analysis including the presence of normal sinus rhythm, Atrial Fibrillation, Bradycardia, Tachycardia and other heart rhythm conditions. Kardia App also allows to collect, analyse, store and share recordings. Quick results in 30 seconds, by simply placing fingers on the two front sensors.

Clinically validated ECG recordings. Unlimited number of ECG tests, anywhere, anytime, without the use of wires, electrodes or gels.

Supplied with CR2016 batteries, quick start guide, Kardia App manual, phone clip and manual in GB, FR, IT.

Compatible with OMRON Connect. Phone clip to attach Kardimobile Stainless steel electrodes to your smartphone.

Technical Specifications

ECG channel: single
 Electrode size: 3x3 cm
 Memory: unlimited
 Compatibility: iOS 10.3.3; Android OS 6.0
 Size: 82x32x3.5 mm
 Weight: 18 g
 Classification: type CF applied part, Class IIa