



Product Code	54101
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2340632
CND	Z12050380
EAN/UPC	4015672114029

## Description

KARDIAMOBILE 6 LEAD ECG - KM-6L-E

KardiaMobile are pocket-sized personal ECGs with finger electrodes that work together with smartphone or tablet through Kardia App.

Visit the following link to view the full list of countries where the product and app are currently available:

<https://alivecor.zendesk.com/hc/en-us/articles/115015719967-Which-Kardia-products-are-available-in-my-Country>

The KardiaMobile System (device + Kardia App) displays ECG rhythms and output of ECG analysis including the presence of normal sinus rhythm, Atrial Fibrillation, Bradycardia, Tachycardia and other heart rhythm conditions. Kardia App also allows to collect, analyse, store and share recordings. Quick results in 30 seconds, by simply placing fingers on the two front sensors. KardiaMobile 6L has an additional electrode for left leg (knee or inside of the ankle) on the back; it delivers six times more accurate data compared to 1L.

Clinically validated ECG recordings. Unlimited number of ECG tests, anywhere, anytime, without the use of wires, electrodes or gels.

Supplied with CR2016 batteries, quick start guide, Kardia App manual and manual in GB, FR, IT.

Compatible with OMRON Connect

## Technical Specifications

ECG channel: six

Electrode size: 3x3 cm

Memory: unlimited

Compatibility: iOS 10.3.3; Android OS 6.0

Size: 89.9x29.9x7.2 mm

Weight: 24 g

Classification: type CF applied part, Class IIa